---

# Chemistry experiments in English

---

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## Measurement of Heat of Reaction

### Objective:

To understand the generation and absorption of heat during chemical reactions by measuring the heat of solution and neutralization of a sodium hydroxide solution.

### Instruments:

You need to prepare a polystyrene cup / graduated cylinder / beaker / a film case / a magnetic stirrer / a stirring bar / a PC / a temperature sensor / a data logger / tweezers

### Reagents:

1.0 mol/L hydrochloric acid, 1.0 mol/L sodium hydroxide solution, and sodium hydroxide (solid)

### Cautions:

- ※ You need to wear safety glasses during the experiment.
- ※ If you touch some reagents with your hand, you need to wash them immediately.

### Procedure:

#### **Experiment 1: Measuring heat of solution of sodium hydroxide**

- I. Measure 100 mL of distilled water using a graduated cylinder and put it into a polystyrene cup with a stirring bar.
- II. Begin measuring the temperature of water inside the cup using a temperature sensor attached to a PC. Measure for 1 minute, and check that the temperature reading is stable.
- III. Put sodium hydroxide (4.0 g, already prepared in a film case) into the cup and dissolve it using a magnetic stirrer. Measure the temperature for 5 minutes.
- IV. Finish measurement and read the change in temperature from the graph
- V. Empty the solution from III into the beaker on the teacher's desk. Rinse the cup, stirring bar and the temperature sensor first with tap water, then with distilled water.

#### **Experiment 2: Measuring heat of neutralization**

- I. Measure 50 mL of 1.0 mol/L hydrochloric acid using a graduated cylinder and put it into a polystyrene cup.
- II. Put a stirring bar into the cup and stir the solution with a magnetic stirrer
- III. Measure the temperature of the solution inside the cup. Measure for 1 minute, and check that the temperature reading is stable.
- IV. Measure 50 mL of 1.0 mol/L sodium hydroxide solution using a graduated cylinder.
- V. Slowly pour the sodium hydroxide solution into the polystyrene cup. Measure the temperature for 5 minutes.
- VI. After finishing measurement, stop the temperature sensor, and read the change in temperature from the graph
- VII. Empty the solution from ⑤ into the sink. Rinse the measuring cylinders, temperature sensor, cup and stirring bar with tap water

### Data Processing

Based on the graphs, analyze how many degrees (K) the temperature rose in Experiment 1 and 2.

Experiment 1 \_\_\_\_\_ K, Experiment 2 \_\_\_\_\_ K

### Experiment 1

- Calculate heat energy generated in this experiment.  
specific heat of solution:  $4.2 \text{ J/(g} \cdot \text{K)}$  density of water:  $1.0 \text{ g/cm}^3$
- Calculate the heat of solution per 1 mol of sodium hydroxide. H=1.0, O=16, Na=23
- Write a thermochemical equation which expresses the heat of solution of sodium hydroxide.

### Experiment 2

- Calculate heat energy generated in this experiment.  
specific heat of solution:  $4.2 \text{ J/(g} \cdot \text{K)}$  density of solution:  $1.0 \text{ g/cm}^3$
- Calculate the heat of neutralization per 1 mol of  $\text{H}^+$  and  $\text{OH}^-$ .
- Write a thermochemical equation which expresses the heat of neutralization.

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